



SPHEROPEEL
— A E S T H E T I C —



What is Spheropeel™?

It's a microdermal
a b r a s o r
by controlled
Radiofrequency

What is it indicated for?

E x f o l i a t i o n
for consequential
c u t a n e o u s
regeneration



Indications:



Peeling



Post-acne scar treatment;



Stretch mark treatment;



Revitalizing treatment: illuminating and hydrating effect;



Treatment of periorcular, glabellar, nasogenien wrinkles with product vehiculation;



Antiaging, sun and age spots.



How does it work™?

Thanks to the combination of controlled power radiofrequency, in synergy with revitalizing and exfoliating substances and characters, depending on the treatment to realize.



How many programs are available?

3 programs:

P E E L I N G
power 3W - frequency 2000 KHZ

V E H I C U L A T I O N
power 4,5W - frequency 1900 KHZ

A N T I A G I N G
power 6W - frequency 900 KHZ



For each program *Spheropeel™* the device is equipped with a different and sterilizable tips (bipolar electrodes):



EAGLE PROBE

For the PEELING program in combination with our Minitech peeling.




SILK PROBE

A big rounded tip for the ANTIAGING program, for rejuvenation in combination with revitalizing substances.



TARGET PROBE

A triangular tip for the VEHICULATION program in combination with revitalizers for the treatment of acne scars, stretch marks and periorcular, glabellar and nasogenien wrinkles.



BEFORE



IMMEDIATELY AFTER



AFTER 1 WEEK

BEFORE



Ultrasound image of the gluteal region with a 7 MHz probe before the application of a revitalizer to vehiculate (fig 1°)

AFTER




Image of the same gluteal region after the vehiculation of the revitalizer

BIBLIOGRAPHY

- 1) Groves, Nancy. "Radiofrequency Device Offers Advantages of its Own" (featuring Brian D. Zelickson, M.D.). *Dermatology Times* Vol. 23, No. 7; July 2002.
- 2) Guttman, Cheryl. "RF Heating Promising for Skin-Tightening, Browlifting" (featuring Richard E. Fitzpatrick, M.D.). *Cosmetic Surgery Times* Vol. 5, No. 7 (02); August 2002.
- 3) Guttman, Cheryl. "Skin Laxity Improvements Seen Vividly in Neck Areas" (featuring Tina S. Alster, M.D.). *Cosmetic Surgery Times*; May 2003.
- 4) Kazanjian, Dodie. "Facelift of the Future" (featuring Roy Geronemus, M.D., Javier Ruiz-Esparza, M.D., Karyn L. Grossman, M.D., M.D., Pat Wexler, M.D.). *Vogue*; May 2003.
- 5) Goldberg DJ, Koch J, Burns J, Geronemus R. "A New Non-Ablative Radio Frequency Device to Tighten Tissue: A Randomized Trial Comparing a Single Treatment to Multiple Treatments".
- 6) Alster TS, Tanzi EL. "Treatment of Prominent Nasolabial Folds and Cheek Laxity with a Nonablative Radiofrequency Device".
- 7) Tanzi EL, Alster TS. "Improvement of Neck Laxity with a Nonablative Radiofrequency Device: A Lifting Experience".
- 8) Fitzpatrick RE, Iyer S. "Collagen and Tissue Tightening: Paradigm for Treatment of the Cheeks and Neck Using Radiofrequency".
- 9) Fitzpatrick RE, Geronemus R, Goldberg D, Klammer M, Kilmer S, Ruiz-Esparza, J. "First Multi-Center Study of a New Non-Ablative Radio Frequency Device to Tighten Facial Tissue".
- 10) Ruiz-Esparza J. "Non-Ablative Radiofrequency Tissue Tightening of Facial Skin: The Medical Face-Lift: A Report of 25 Patients".
- 11) Alster TS. "Radiofrequency Treatment: A Non-Surgical Approach to Skin Tightening and Other Applications".
- 12) Ruiz-Esparza J. "The Non-Surgical Face-Lift Via Non-Ablative Radiofrequency".
- 13) Nahm W. "Radiofrequency applied to Half of Face Ipsilaterally Produces Elevation of The Brow, Changes in the Angle of the Eyebrow, and Improvement in Appearance of Melolabial Fold and Jawl".
- 14) Fitzpatrick RE, Geronemus R, Goldberg D, Klammer M, Kilmer S, Ruiz-Esparza, J. "Multicenter Study of Noninvasive Radiofrequency for Periorbital Tissue Tightening". *Lasers in Surgery and Medicine* Vol. 33, No. 4; 2003.
- 15) Iyer S, Suthamjaryiya K, Fitzpatrick RE. "Using a Radiofrequency Energy Device to Treat the Lower Face: A Treatment Paradigm for a Nonsurgical Facelift". *Cosmetic Dermatology* Vol. 16, No. 2; February 2003.
- 16) Ruiz-Esparza J, Gomez JB. "The Medical Face Lift: A Noninvasive, Nonsurgical Approach to Tissue Tightening in Facial Skin Using Nonablative Radiofrequency". *Dermatol Surg* Vol. 29, No. 4, p 325-332; April 2003.
- 17) Zelickson BD. "Histological and Ultrastructural Evaluation of the Effects of a Radiofrequency-Based Nonablative Dermal Remodeling Device". *Archives of Dermatology* Vol. 140, p 204-209; February 2004.
- 18) Alster T.A., Tanzi E. "Improvement of Neck and Cheek Laxity With a Non-Ablative Radiofrequency Device: A Lifting Experience". *Dermatol Surg* Vol 30, p 503-507; April 2004. Book Chapters

PROMOITALIA

WELLNESS RESEARCH

ITALIAN BRANCHES

Milan 20124
Via San Gregorio, 44 - ITALY
Ph. +39 02 670 77 226
Fax +39 02 669 89 000

ASIAN BRANCH

Bangkok 10500
12th Floor, Bubhagit Building
20 North Sathorn Road - THAILAND

EXPORT OFFICE

Naples 80078
Via Antiniana, 59 - ITALY
Ph. +39 081 55 85 118
Fax +39 081 22 92 267

Rome 00195
Via F. De Sanctis, 15 (Prati) - ITALY
Ph. +39 06 375 11 826 - Showroom
Ph. +39 06 972 743 52
Fax. +39 06 972 743 10

PROMOITALIA MIDDLE EAST FZ LLC

E Block ,Ground Floor,
Building No. 64
Dubai Healthcare City, Dubai U.A.E

www.webpromoitalia.com

beauty.instruments@webpromoitalia.com